



Govt. of West Bengal
Office of the Principal
College of Medicine & Sagore Dutta Hospital
Kamarhati, Kolkata-700058

Ph no: 033-2583-4277

Memo No: CMSDH/Pr./0072/18

Date: 15-01-2018

Quotation Notice

The quotation is invited from the reputed organization for the supply of following equipments at **Department of Pathology** under College of Medicine & Sagore Dutta Hospital.

<u>Sl.</u>	<u>Name of the equipments</u>	<u>Sepecification</u>	<u>Quantity</u>
01.	Digital tissue floatation bath	As per annexure-I	As required
02.	Bone saw	As per annexure-II	As required

The rate should be inclusive of all taxes/GST.


The following documents should be submitted along with quotation letter.

01. Xerox copy of Trade license.
02. Xerox copy of VAT & GST registration.
03. Xerox copy PAN card.
04. Last year IT Return (2015-16)/(2016-17).

The interested organization may communicate for specification.

The quotation must reach this office on or before **29.01.2018**, with super scribed on the top of the envelope. Quotation will be opened on **30.01.2018** at 1 P.M. at the office of the undersigned.

The undersigned reserves the right to accept or reject any quotation without assigning any reason.


Principal


College of Medicine & Sagore Dutta Hospital
Kamarhati, Kolkata-700058

Memo No: - CMSDH/Pr./0072/18

Date: 15-01-2018

Copy forwarded for necessary information:-

1. MSVP, College of Medicine & Sagore Dutta Hospital, Kamarhati, Kolkata-700 058.
2. Accounts Officer, College of Medicine & Sagore Dutta Hospital, Kolkata-700 058.
3. Notice Board, College of Medicine & Sagore Dutta Hospital, Kamarhati, Kolkata-700 058.
4. HOD, Pathology, College of Medicine & Sagore Dutta Hospital, Kamarhati, Kol-700 058.


Principal

College of Medicine & Sagore Dutta Hospital
Kamarhati, Kolkata-700058

Annexure-I

Digital tissue floatation bath

- Inner tank (single piece) made of black anodized aluminium.
- Housing made of powder coated aluminium.
- Front panel should consist of temperature regulator knob, mains switch and pilot lamp.
- Inside black and outer surface finished in powder paint, with wide enough rim for drying the wet slides.
- Thermal loss should be prevented by filling mineral wool (insulation) in space between the outer body & inner tank.
- Heating elements should be positioned in the floor and side walls and a natural mixing of the medium arises, providing better temperature distribution in the bath.
- Temperature controlled by imported capillary type thermostat, from ambient to 70°C with an accuracy of +2°C.
- Low 110mm depth of the bath to allow comfortable and safe working.
- Should be supplied with Dust Guard Lid made of black anodized aluminium.
- Operating temperature range: +15 °C to +40 °C
- Temperature control range: Room temperature to 75 °C